



1

HYPERLAB METHODS MANUAL

April 2020 Revision



STEROGLASS S.r.l.
apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.
Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia
Tel. 075.609091 (r.a.) Fax 075.6090950
e-mail: info@steroglass.it ; steroglass.amm@pec.collabra.it ; internet: www.steroglass.it
R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Preliminary Instructions

In some reaction kits it is possible to find the BLANK reagent; this product is almost never called for in the performing of the Hyperlab methods applied to robotic instrumentation. BLANK reagent is found in certain reaction kits because they can also be used with manual or automatic flow cell instruments where the determination of the sample blank must necessarily be carried out separate from that of the reagent blank.

Once the daily analysis cycle has been completed, the bottles containing the reagents must be capped and placed in a refrigerator (4° to 8°C) together with the reagent holder. Leaving the bottles uncapped on the cooling plate in operation brings about a greater adsorption of oxygen by the reagents, which changes their titer, and consequently may compromise the results of subsequent analyses.

For correct use of Hyperlab, the water supply tank must be filled with double-distilled or distilled water, adding 1 ml of *System Solution* for each liter of water used. For example, when using 15 l of water, add 15 ml of system solution. The washing tank, provided exclusively for Hyperlab Plus, must instead be filled with a 0.1 N sodium hydroxide (NaOH) solution using double-distilled or distilled water. For example, to make 2 l of this solution, weigh 4 g of pure NaOH (drops or flakes) and dissolve it in 2 l of water.

For Hyperlab Plus and Basic it is recommended that the cuvettes cleaning procedure be carried out once a week from the *Maintenance Procedures* menu using the alkaline solution in position 1 (1-part bleach and 3 parts water) and the acid solution in position 2 (hydrochloric acid 0.1 N). When performing the methods that are most "dirtying," such as polyphenols, catechins, sulfur dioxide, chlorides and calcium, the cleaning procedure should be carried out once per day. For Hyperlab Smart instead it is recommended that the sampling needle be washed monthly.

The kit components are stable until the expiration date indicated on the package. They must be stored away from light and oxygen at 4° to 8°C (in the refrigerator).

Whenever possible, it is recommended that analyses be performed with kits and Auto methods that have a ratio of 4:1 between reagent 1 and reagent 2 instead of 40:1 as for traditional methods, and there is no blank reagent. These methods are more efficient and less expensive as regards managing reagents.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@steroglass.it ; steroglass.amm@pec.collabra.it ; internet: www.steroglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



CONTENTS

ACETIC ACID AUTO	4
ACETIC ACID	4
ACETALDEHYDE	5
ANTHOCYANINS	6
RAN.....	6
ASCORBIC ACID	6
CALCIUM	7
CATECHINS	7
CITRIC ACID	8
CHLORIDES	8
COLOR	9
COPPER SELF BLANK	9
D-LACTIC ACID	10
D-MALIC ACID	10
OD 420 nm.....	11
OD 520 nm.....	11
ETHANOL.....	12
FRUCTOSE	12
BIVALENT IRON	13
IONIC IRON	13
GLUCOSE AUTO.....	14
GLUCOSE + FRUCTOSE	14
GLUCOSE FRUCTOSE 0-7 g/l AUTO	15
GLUCOSE FRUCTOSE 10-50 g/l AUTO	15
GLUCOSE FRUCTOSE 50-300 g/l AUTO	16
GLUCOSE FRUCTOSE 0-7 g/l	16
GLUCOSE FRUCTOSE SUCROSE	17
GLYCERIN	17
UV GLYCEROL	18
GLUCONIC ACID	18
POTASSIUM	19
AUTO L-LACTIC ACID	19
L-LACTIC ACID	20
AUTO L-MALIC ACID	20
L-MALIC ACID	21
MAGNESIUM.....	21
AMINO NITROGEN	22
AMMONIACAL NITROGEN	22
PYRUVIC ACID	23
WHITE WINE POLYPHENOLS	23
TOTAL POLYPHENOLS	24
PROTEINS.....	24
SUCROSE	25
FREE SULFUR DIOXIDE	25
TOTAL SULFUR DIOXIDE	26
TARTARIC ACID	27



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



The Hyperlab methods for the Plus, Basic and Smart models and the characteristics of the reagent kits for each individual analysis method are described below according to the following scheme:

METHOD CODE – METHOD NAME – REAGENT KIT CODE NO. – REAGENT KIT VOLUME

The description of each individual method includes an introductory part that gives the following information about the analyte: origin, oenological function, expected concentration and legal terms (if applicable). This is followed by a table regarding the use of the reagents found in the reaction kit and instructions for preparing the reagents when necessary.

The reaction principles and the characteristics of the method are then indicated.

Lastly, if there is the reference *Notes*, recommendations and tips are provided for optimizing the analysis.

ACETAU – ACETIC ACID AUTO – SQPE068205 – 125 ml

It's a product of oxidation of ethanol and represents about 95% of volatile acidity in wine. The volatile acidity provides information on the health of a wine, whether legacy of disease, winemaking wrong or faulty storage.

The optimum concentration varies from 0.2 to 0.7 g / l (e.g. legal limit for many standards is about 1 g / l, dependent on the degree of alcohol).

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
ACET AUTO R1A (ACA1A)	1 x 50	large
ACET AUTO R1B (ACA1B)	1 x 50	large
ACET AUTO R2 (ACAR2)	1 x 25	small

Method principle and characteristics

Acetyl-CoA-synthetase converts acetic acid into acetyl-CoA, which reacts with oxaloacetate in the presence of citrate synthase to form citrate. The oxaloacetate necessary for the reaction is produced through the activity of the enzyme malate dehydrogenase on malate and on NAD, therefore NADH is generated. The increase in absorbance due to the formation of NADH, measured at 340 nm, is proportional to the amount of acetic acid in the sample.

Linearity with total V / sample V = 75: up to 1.2 g/l

Precision: CV% <1.5

Calibration and control: multiparameter standard SQPE053234 acetic acid 1g/l; multilevel standard SQPE081638 acetic acid 1 g/l (level 1), 0.8 g/l (level 2), 0.6 g/l (level 3), 0.4 g/l (level 4), 0.2 g/l (level 5).

ACETIC – ACETIC ACID – SQPE059575 – 5 x 20 ml

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
ACET R1A (ACE1A)	5 x 10	small
ACET R1B (ACE1B)	5 x 10	small
ACET R2 (ACER2)	1 x 2.5	cup
ACET BL	2 x 50	do not use

Method principle and characteristics

Acetyl-CoA-synthetase converts acetic acid into acetyl-CoA, which reacts with oxaloacetate in the presence of citrate synthase to form citrate. The oxaloacetate necessary for the reaction is produced through the activity of the enzyme



STEROGLOSS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



malate dehydrogenase on malate and on NAD, therefore NADH is generated. The increase in absorbance due to the formation of NADH, measured at 340 nm, is proportional to the amount of acetic acid present in the sample.

Linearity with total V / sample V = 75: up to 1.2 g/l

Precision: CV% <1.5 in white wine; CV% <1.75 in red wine

Calibration and control: multiparameter standard SQPE053234 acetic acid 1g/l; multilevel standard SQPE081638 acetic acid 1 g/l (level 1), 0.8 g/l (level 2), 0.6 g/l (level 3), 0.4 g/l (level 4), 0.2 g/l (level 5).

ALDEID – ACETALDEHYDE - SQPE059576 – 5 x 20 ml

Acetaldehyde is the product of intermediate oxidation of ethyl alcohol, before this becomes acetic acid. Being index of initial oxidation, it is good to determine the acetaldehyde throughout the processing phase and also in the refining phase, as well as before bottling. The determination of acetic aldehyde is also very important during the practice of micro-oxygenation (its concentration remains constant)

5

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
ALD R1A (ALD1A)	5 x 20	small
ALD R1B (ALD1B)	1 x 2.5	cup
ALD R2 (ALDR2)	1 x 2.5	cup
ALD BL	2 x 50	do not use

Homogenize the ALD R1B suspension by shaking manually before use.

Method principle and characteristics

Acetaldehyde is quantitatively oxidized by NAD to acetic acid in the presence of the enzyme aldehyde dehydrogenase (ALDH).

The amount of NAD reduced to NADH during the reaction, measured at 340 nm, is directly proportional to the amount of acetaldehyde present in the sample.

Linearity with total V / sample V = 20: up to 150 mg/l

Precision: CV% < 1.5 on decolorized sample

Calibration and control: standard SQPE053229 acetaldehyde 100 mg/l.

Notes: With the repetition of the sample/standard over time, it is advisable to renew the sample, given the high volatility of the acetaldehyde. In addition, to improve the accuracy of the data, red wines can be decolorized, especially if the total polyphenols are over 2500 mg/l, using activated carbon or PVPP 1% weight/volume.

ANTOCI – ANTHOCYANINS - SQPE054971 – 40 x 50 ml

They are part of the phenolic compounds, the main red and blue dyes and extremely important to the organoleptic characteristics of a wine. They are determined during the maceration of red grapes, during and after fermentation. Together with the dosage of total polyphenols help to determine whether a product deserves to be "refined".

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
ANTO R1 (ANTOC)	4 X 50	large



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Method principle and characteristics

Anthocyanins are ionized in an acidic environment by a controlled ionic strength buffer, which also contains solubilizers for eliminating any turbid protein formations, and measured at 520 nm. The anthocyanin fraction involved in polymerization with other polyphenolic substances is not detected.

Linearity with total V / sample V = 10: up to 600 mg/l

Precision: CV% < 1.5

Result unit of measurement: mg/l

Notes: Must that is excessively turbid needs to be centrifuged or filtered before the analysis.

RAN – READILY ASSIMILABLE NITROGEN

R.A.N. is the sum of AMMONIACAL NITROGEN and ALPHA-AMINO NITROGEN. The available nitrogen (RAN) is also a factor of great importance to determine whether there are sufficient nutrients to the growth of yeasts responsible for fermentation. They are really important in relation to protein stability and evaluation of the "authenticity" of the wine: ammonia nitrogen may be added to mask the lack of amino acids derived from grapes.

This is a compound method and is derived from the sum of the concentration of amino nitrogen and the concentration of ammoniacal nitrogen. The characterization of the reagents and of these two methods is indicated below. In Hyperlab code terms: $NH_2 + NH_3 = RAN$ (mg/l).

ASCORB – ASCORBIC ACID - SQPE072166 – 115 ml

The L-ascorbic acid, vitamin C, is a naturally occurring organic compound with antioxidant properties widely used on food & beverage industry to prevent oxidation of aromatic fractions and coloring of various beverages such as wines, fruit juices, etc

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
L-ASC R1 (ASCR1)	2 X 50	large
L-ASC R2 (ASCR2)	1 x 10	small
L-ASC R3 (ASCR3)	1 x 5	cup

These 3 reagents are the same used in the ASCR and ASCB coded methods that start automatically launching the ASCORB method.

Method principle and characteristics

The concentration of ascorbic acid in the sample is provided by the difference between the total quantity of reducing agents in the sample (ASCR) and a blank sample in which the ascorbic acid (ASCB) is removed. Therefore $ASCR - ASCB = ASCORB$, or the concentration of ascorbic acid. These determinations are performed automatically by launching exclusively the ASCORB method.

At an enzymatic level, all the reducing agents present in the sample reduce MTT (Thiazolyl Blue Terazolium Bromide) in the presence of PMS (Phenazine methylsulfate) to form formazan. The removal of ascorbic acid, upon which the calculation is based, is instead carried out by ascorbate oxidase. The difference in absorbance obtained, measured at 578 nm, is proportional to the amount of L-ascorbic acid present in the sample.

Linearity with total V / sample V = 50: up to 300 mg/l

Precision: CV% < 1.0

Calibration and control: standard preparation at 250 mg/l, dissolve 125 mg/l of ascorbic in a 500 ml flask and then bring to volume.

Notes: Add 0.5 ml of propionaldehyde with a concentration of 10 g/l in 2 ml of sample to limit interference from SO_2 . Wait 5 minutes before starting the analysis. The result obtained must be multiplied by 1.25 ($0.5 \text{ ml} + 2 \text{ ml} = 2.5 \text{ ml}$; $2.5 \times 0.5 = 1.25$). Do not leave the sample to be analyzed too long on the sample holder, given the high reactivity of ascorbic acid.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Ca²⁺ – CALCIUM – SQPE054972 – 2 X 100 ml

In juices tends to precipitate as for example during alcoholic fermentation in wine.

It's important to determine the concentration (which should not exceed about 80 g/l) to prevent precipitation of neutral calcium tartrate.

The colorimetric method responds very well to the need for this determination.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
CA R1 (CALCR)	2 x 100	small
CA BL (CALCB)	2 x 100	small

7

Method principle and characteristics

In a sub-acidic environment, calcium reacts with Arsenazo III, a highly selective chromogen for this element, producing a colored complex of intensity proportional to the amount of calcium present in the sample. Optical density is measured at 650 nm.

Linearity with total V / sample V = 60: up to 150 mg/l

Precision: CV% < 2.5

Calibration and control: standard SQPE053231 calcium 100 mg/l

Notes: Centrifuge or filter excessively turbid must before the analysis. Given the high sensitivity of this method, it is best to avoid any calcium contamination that may derive from: water that is not perfectly distilled and/or demineralized or from glassware with calcium residues. Furthermore, before performing the analysis, it would be best to carry out the "cuvette cleaning" procedure with two bottles containing 0.1N hydrochloric acid, in addition to performing this method on its own.

CATEC – CATECHINS – SQPE054972 – 5 x 20 ml

Catechins are part of the family of polymers condensed tannins present in grapes and wine. The analysis is carried out during the period of maceration, on must, before and after fermentation on the finished wine.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
CATE R1 (CATRM)	5 x 15	small (3 parts)
CATE R2 (CATRM)	5 x 5	small (1 part)
CATE BL	2 x 50	do not use

Mix three parts of R1 with one part of R2 to obtain a single CATRM working reagent (small bottle). It is possible to add one bottle of R2 to one of R1 (1:3), making sure to store any remaining amount at 2-8°C away from light.

Method principle and characteristics

In a strongly acidic environment and in a non-aqueous solvent, catechins react with cinnamaldehyde to form a chromophore with maximum absorption at 644 nm. Consequently, the measurement is done at 650 nm.

Linearity with total V / sample V = 6: up to 50 mg/l

Precision: CV% < 1.5

Calibration and control: standard SQPE082087 catechin 500 mg/l (dilute 10 times to obtain the standard at 100 mg/l). Once diluted it will keep for one week in the refrigerator.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



CITRC – CITRIC ACID – SQPE0076313 – 5 x 10 ml

Organic acid naturally present in several juices. Contributes to the formation of total acidity and has the property to involve ferric iron in a soluble complex anion. The analysis is performed to determine the concentration present and evaluate possible additions to reach the desired level of total acidity.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
CITR R1A (CIT1A)	5 x 10	small
CITR R1B (CIT1B)	1 x 1.25	cup
CITR R2 (CITR2)	5 x 0.25	cup
CITR BL	1 x 50	do not use

The R2 reagent (lyophilized) must be reconstituted by pouring 0.5 ml of distilled water into the cup. It remains stable for about 2 days. Shake the R1B reagent well before use.

Method principle and characteristics

In the presence of the enzyme citrate lyase, citric acid is converted to oxaloacetate, which is then reduced to malate by malate dehydrogenase, in the presence of NADH. Due to the possible decomposition of the oxaloacetate in reaction conditions, any pyruvate formed is converted to lactate by lactate dehydrogenase, in the presence of NADH. The decrease in absorbance due to the consumption of NADH, measured at 340 nm, is proportional to the amount of citric acid present in the sample.

Linearity with total V / sample V = 75: up to 0.8 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 citric acid 1g/l; multilevel standard SQPE081638 citric acid 0.8 g/l (level 1), 0.6 g/l (level 2), 0.4 g/l (level 3), 0.3 g/l (level 4), 0.1 g/l (level 5).

Notes: In the case of strongly colored red wine (total polyphenols > 2500 mg / l) it is advisable to decolorize the sample with activated carbon or 1% PVPP weight/volume.

CL – CHLORIDES – SQPE055024 – 2 x 100 ml

Chloride are present in any food stuffs in varying amounts.

The analysis on the finished product is carried out to verify the content according to the law limits.

Analysis is usually carried out by titration while colorimetric method, which is also very accurate, is now in routine use of many laboratories.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
CL R1 (CLORR)	2 x 100	small
CL BL (CLORB)	2 x 100	small

Store the reagents at room temperature (15-20°C) and away from light.

Method principle and characteristics

Chloride ions react with mercury thiocyanate, developing a proportional amount of thiocyanate ions. The thiocyanate ions react with the iron (III) ions present in solution to form a red-colored complex characterized by an absorbance peak at 480 nm. The optical density is thus measured at 492 nm.

Linearity with total V / sample V = 60: up to 1.5 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE056370 chlorides 1 g/l NaCl



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



COLORE – COLOR – SQPE054875 – 4 x 100 ml

Really important parameter to qualify a wine and its commercial value.

The reading of the samples at 420 and 520 nm determines the color points (intensity) This method also allows even reading at 620 nm (dark tones).

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
COL R1 (COL1)	4 x 100	large

Method principle and characteristics

The sample in question undergoes various programmed dilutions by means of a buffer stabilized at a known pH and with high ionic strength, thus avoiding the alteration of the ratio between the substances that contribute to the formation of the color. Following this, a mathematical algorithm is used to compare the results to the sample as is, therefore without considering the dilutions.

Notes: to observe the complete results, click on the green bar of the test of the sample being examined without dilution type C (Completed), which gives: the color intensity **IC**, the modified color intensity **ICM**, the hue **TON** and the pie chart with the different % contributions of the three absorbance values that contribute to the color formation.

Linearity up to 40 OD

Precision: CV% < 1.5 in red wine

Cu2+ – COPPER SELF BLANK – SQPE075544 – 50 ml

Juices may contains copper both from fruit (minimally) and from cupric treatments. The copper should be determined after fermentation to assess whether or not to proceed to the stage of "de-metallization", as there is a legal limit expected.

In addition, monitoring of the copper concentration is considered necessary in the pre-bottling stage, since it is subject to a legal limit of 1 mg/l.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
CU R1	1 x 21	small
CU R2	1 x 4	cup
CU BL	1 x 25	use as R1

The reagent kit must be stored outside the refrigerator at room temperature (15±25°C). If it is at a lower temperature, the R1 and/or BL reagent must be warmed, for example in a thermostated bath at a temperature of about 30°C.

Method principle and characteristics

The 3,5-Di-Br-PAESA chromogen reacts with copper 2+ in the presence of sodium dodecyl sulfate to form a blue-violet complex, the absorbance of which is measured at 578 nm. The amount of chromophore formed is proportional to the amount of copper present in the sample. The specificity of the reaction is ensured by a specific environment, especially as regards the pH value.

Linearity with total V / sample V = 3: up to 1.2 mg/l

Precision: CV% < 2

Calibration and control: standard SQPE056386 copper 2 mg/l.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

D-LATT –D-LACTIC ACID – SQPE059194 – 5 x 20 ml

Thus the presence of D-lactic acid (D+ dextrorotatory) indicates an unpleasant bacterial development starting from the residual sugars which always goes along with the production of acetic acid. In wine this fault is called mannitic fermentation.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
D-LAT R1 (DLAT1)	5 x 20	small
D-LAT R2 (DLAT2)	1 x 2.5	cup
D-LAT BL	2 x 50	do not use

Method principle and characteristics

D-lactate is converted into pyruvate by the enzyme D-lactate dehydrogenase in the presence of NAD. Given that the reaction equilibrium is clearly in favor of the D-lactate, the pyruvate that is formed is progressively removed from the reaction environment by the activity of the enzyme D-glutamate pyruvate transaminase, which transforms it into D-alanine.

The NADH that is formed from the oxidation of D-lactate is in a stoichiometric concentration with the analyte D-lactic acid. Thus the absorbance is measured at a wavelength of 340 nm.

Linearity with total V / sample V = 75: up to 1 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 D-lactic acid 3 g/l; multilevel standard SQPE081638 D-lactic acid 3.0 g/l (level 1), 2.4 g/l (level 2), 1.8 g/l (level 3), 1.2 g/l (level 4), 0.6 g/l (level 5).

D-MAL – D-MALIC ACID – SQPE067017 – 5 x 25 ml

The D(+) dextrorotatory isomer of malic acid may be present in traces; significant quantities can arise from food adulteration. The addition of D-malic acid for correcting the acidity of wines is currently allowed by law.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
D-MAL R1A (DML1A)	5 x 20	small
D-MAL R1B (DML1B)	1 x 25	small
D-MAL R2 (DMLR2)	1 x 2.5	cup
D-MAL BL	2 x 62.5	do not use

Shake R2 well before using.

Method principle and characteristics

The enzyme D-malate dehydrogenase catalyzes the selective oxidation of D-malic acid in the presence of NAD. The isomer L-malic acid, naturally present in solution, does not undergo reaction. The increase in absorbance due to the formation of NADH, measured at 340 nm, is proportional to the amount of D-malic acid in the sample.

Linearity with total V / sample V = 15: up to 1 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 D-malic acid 5 g/l; multilevel standard SQPE081638 D-malic acid 5.0 g/l (level 1), 4.0 g/l (level 2), 3.0 g/l (level 3), 2.0 g/l (level 4), 1.0 g/l (level 5).

DO420 – OD 420 nm – SQPE054875 – 4 x 100 ml

The color of wine is one of the main parameters for sensory evaluation. It also reflects some characteristics of the product such as the grape variety, the vintage, and the evolutionary state. In oenology it is common practice to observe absorbance values at 420 nm, 520 nm and 620 nm, since the sum 420 nm + 520 nm gives the coloring intensity of the product, and further adding the absorbance value at 620 nm one obtains the modified coloring intensity, while the 420 nm /520 nm ratio provides the hue of the product. The direct reading on the sample of the optical density at 420 nm can be useful in the case of unwanted oxidation processes in which, at the chromatic level, the contribution given by the wavelength at 420 nm increases.

This method does not provide for the use of reagents.

Method principle and characteristics

The sample being examined is read spectrophotometrically at 420 nm.

Linearity up to 3 OD

Notes: direct reading method suitable for white and rosé wines.

DO520 – DO 520 nm – SQPE054875 – 4 x 100 ml

The color of wine is one of the main parameters for sensory evaluation. It also reflects some characteristics of the product such as the grape variety, the vintage, and the evolutionary state. In oenology it is common practice to observe absorbance values at 420 nm, 520 nm and 620 nm, since the sum 420 nm + 520 nm gives the coloring intensity of the product, and further adding the absorbance value at 620 nm one obtains the modified coloring intensity, while the 420 nm /520 nm ratio provides the hue of the product. The direct reading on the sample of the optical density at 520 nm can be useful for evaluating the red component of the color of the wines, which tends to change according to the aging techniques used.

Method principle and characteristics

The sample being examined is read spectrophotometrically at 520 nm.

Linearity up to 3OD

Notes: direct reading method suitable for white and rosé wines.

ETANOL – ETHANOL – SQPE078690 – 50 ml

Ethyl alcohol is produced through fermentation processes carried out by yeasts that convert sugars into carbon dioxide and ethanol, its concentration increases as the fruit ripens and has to be checked on finish juice as per legal limit (0,5% Vol on fruit juices).

This kit can be used for low ethanol concentration, lower than 2,5%Vol, both in juices than at initial alcoholic fermentation in wine.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
ETAN R1 (ETAN1)	1 x 25	large
ETAN R2 (ETAN2)	1 x 25	small

Method principle and characteristics

Ethanol is oxidized to acetaldehyde by alcohol dehydrogenase in the presence of NAD, resulting in the formation of NADH. The NADH that is formed is in stoichiometric concentration with ethanol, therefore the increase in NADH is measured spectrophotometrically at 340 nm.

Linearity with total V / sample V = 100: up to 2,5% vol.

Precision: CV% < 1.5 up to 2,5% vol.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

F – FRUCTOSE – SQPE063151 – 5 x 20 ml

Fructose is one of the simple sugars of which there is the most in fruit. It differs from glucose in that it has a ketonic function (less reactive) and a greater sweetening power as regards the sensory analysis

A few g/l of fructose in wine can represent the energy source for any bacterial development; in finished products it gives sweetness and roundness. It is more appropriate to discriminate simple sugars individually, and the enzymatic method is identified as official for these determinations.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GL+FR R1 (G+FR1)	5 x 20	small
GL+FR R2 (G+FR)	1 x 2.5	cup
GL+FR R3 (G+FR3)	1 x 2.5	cup
GL+FR BL	2 x 50	do not use

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of fructose present in the sample.

Linearity with total V / sample V = 100: up to 2 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE063151 glucose-fructose-sucrose 20 g/l, with fructose 5 g/l, glucose 5 g/l and sucrose 10 g/l.

Notes: Method indicated for medium-low fructose concentrations. Given that in the presence of the enzyme hexokinase the glucose in the sample is also converted into glucose-6-phosphate, this method is reliable exclusively for samples in which fructose prevails over glucose, which must be decidedly less (fructose/glucose > 2).

Furthermore, the fructose concentration is obtained from the difference between the determination of glucose-fructose minus glucose.

Fe²⁺ - BIVALENT IRON – SQPE062468 – 5 x 20 ml

Fruit juice and wine contains iron: from grapes, from the ground, from manipulation, from transport etc.

The determination is performed to determine whether to proceed to the stage of "de-metallization". Colorimetric method is becoming "routine" as allows precise determination at low concentrations. Measuring Iron ion 2 and 3 allows even to evaluate redox potential in wine.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
FE R1 (FERR1)	5 x 20	small
FE R2	1 X 20	do not use
FE R3 (FERR3)	1 x 2.5	cup
FE BL	2 x 50	do not use

Method principle and characteristics

The Fe (II) ions react with the Ferene-S to give a blue colored complex, whose absorbance at 578 nm is directly proportional to the concentration of iron in the sample. Thiourea is also added to eliminate copper interference.

Linearity with total V / sample V = 14: up to 40 mg/l



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

Precision: CV% < 2

Calibration and control: standard SQPE053230 iron 20 mg/l.

FeION – IONIC IRON – SQPE062468 – 5 x 20 ml

Bivalent iron (II) is a fundamental cofactor of the chemical oxidation, while ionic iron (III), which can derive from the oxidation of bivalent iron as a function of the system's redox potential, can give rise to precipitation. To counter these phenomena, in addition to the use of antioxidants and protective colloids, it is possible to make use of the chelating action of citric acid against iron (III).

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
FE R1 (FERR1)	5 x 20	small
FE R2 (FERR2)	1 X 20	small
FE R3 (FERR3)	1 x 2.5	cup
FE BL	2 x 50	do not use

Method principle and characteristics

The iron is ionized in an acid medium and then reduced by ascorbic acid to Fe (II). The Fe (II) ions react with Ferene-S to give a blue colored complex, whose absorbance at 578 nm is directly proportional to the concentration of iron in the sample. Thiourea is also added to eliminate copper interference.

Linearity with total V / sample V = 14: up to 40 mg/l

Precision: CV% < 2

Calibration and control: multiparameter standard SQPE053230 iron 20 mg/l.

G – GLUCOSE AUTO – SQPE079100 – 50 ml

The sugar concentration is determined to establish the technological maturity of the fruits, to follow the fermentation process and to determine the residual sugars.

Glucose is the simple sugar most present in fruit and is the elective substrate of the energy metabolism of yeasts and many other microorganisms. It differs from fructose in that it has an aldehyde function (more reactive) and a lower sweetening power as regards the sensory analysis. It also represents the aglycone for many volatile aromatic compounds and for anthocyanins.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GL AUTO R1 (GLAU1)	1 x 40	large
GL AUTO R2 (GLAU2)	1 x 10	small

Method principle and characteristics

The glucose is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose present in the sample.

Linearity with total V / sample V = 100: up to 8 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE063151 glucose-fructose-sucrose 20 g/l with glucose 5 g/l, fructose 5 g/l and sucrose 10 g/l.

Notes: Method indicated for medium-low glucose concentrations.



G+F – GLUCOSE + FRUCTOSE – SQPE063019 – 5 x 20 ml

Glucose, fructose are linked the sweetness and energy power of a beverage and, as it ferments, to the alcoholic strength of a wine.

In wine the determination of glucose / fructose is carried out mainly to follow the process of fermentation of wine and to determine the residue at the end of fermentation. Enzymatic way is now official method for many organizations as OIV and IFU. Moreover, Fehling method is quite impractical and also not determines pentose sugars fermentable as sucrose.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GL+FR R1 (G+FR1)	5 x 20	small
GL+FR R2 (G+FR2)	1 x 2.5	cup
GL+FR R3 (G+FR3)	1 x 2.5	cup
GL FR BL	2 x 50	do not use

14

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter, together with the glucose present in the sample, is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose/fructose present in the sample.

Linearity with total V / sample V = 100: up to 8 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE053233 glucose-fructose 20 g/l.

Notes: Method indicated for medium-low glucose-fructose concentrations.

GF-7 – GLUCOSE FRUCTOSE AUTO 0-7 g/l – SQPE068207 – 125 ml

The residual sugars, on the order of a few g/l, give sweetness to the product but at the same time they constitute the energy source for any unwanted bacterial development. This method allow G/F measurement from 0 to 7g/l

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GF AUTO R1 (GFAU1)	1 x 40	large
GF AUTO R2 (GFAU2)	1 x 10	small

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter, together with the glucose present in the sample, is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose/fructose present in the sample.

Linearity with total V / sample V = 100: up to 8 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE053233 glucose-fructose 20 g/l.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Notes: Method indicated for low glucose-fructose concentrations.

GF-50 – GLUCOSE FRUCTOSE AUTO 7-50 g/l – SQPE068207 – 125 ml

The residual sugars, on the order of a few g/l, give sweetness to the product but at the same time they constitute the energy source for any unwanted bacterial development. This method allow G/F measurement from 7 to 50 g/l

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GL AUTO R1 (GFAU1)	1 x 40	large
GL AUTO R2 (GFAU2)	1 x 10	small

15

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter, together with the glucose present in the sample, is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose/fructose present in the sample.

Linearity with total V / sample V = 75: up to 50 g/l

Precision: CV% < 1 in white wine; CV% < 1.5 in red wine

Calibration and control: standard SQPE053233 glucose-fructose 20 g/l; multiparameter standard SQPE067024 glucose/fructose 50 g/l.

Notes: Method indicated for medium glucose-fructose concentrations.

GF-300 – GLUCOSE FRUCTOSE AUTO 50-300 g/l – SQPE068207 – 125 ml

Glucose, fructose are linked the sweetness and energy power of a beverage and, as it ferments, to the alcoholic strength of a wine. This method allow G/F measurement from 50 to 300 g/l

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GF AUTO R1 (GFAU1)	2 x 50	large
GF AUTO R2 (GFAU2)	1 x 25	small

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter, together with the glucose present in the sample, is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose/fructose present in the sample.

Linearity with total V / sample V = 50: up to 300 g/l

Precision: CV% < 2

Calibration and control: multiparameter standard SQPE067024 with six solutions with different glucose-fructose concentrations: 50, 100, 150, 200, 250 and 300 g/l.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

**GFDRY – GLUCOSE FRUCTOSE 0-7 g/l – SQPE053688 – 5 x 20 ml**

The residual sugars, on the order of a few g/l, give sweetness to the product but at the same time they constitute the energy source for any unwanted bacterial development. This method allow G/F measurement from 0 to 7g/l

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GL-FR R1 (GFR1)	5 x 20	small
GL-FR R2 (GFR2)	1 x 2.5	cup
GL-FR BL	2 x 50	do not use

Method principle and characteristics

In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer. The latter, together with the glucose present in the sample, is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose/fructose present in the sample.

Linearity with total V / sample V = 100: up to 8 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE053233 glucose-fructose 20 g/l.

Notes: Method indicated for low glucose-fructose concentrations.

GFS – GLUCOSE FRUCTOSE SUCROSE – SQPE068207– 125 ml

Glucose, fructose and sucrose are linked the sweetness and energy power of a beverage and, as it ferments, to the alcoholic strength of a wine.

In wine the determination of glucose / fructose is carried out mainly to follow the process of fermentation of wine and to determine the residue at the end of fermentation. Enzymatic way is now official method for many organizations as OIV and IFU. Moreover, Fehling method is quite impractical and also not determines pentose sugars fermentable as sucrose.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GF AUTO R1 (GFAU1)	1 x 40	large
GF AUTO R2 (GFAU2)	1 x 10	small
SACC	10 ml	small

The reagent SACC is sold individually with its own packaging, identified by the code SQPE063020 – SUCROSE – 10ml.

Method principle and characteristics

Sucrose is hydrolyzed by the enzyme invertase into D-glucose and D-fructose. In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer, which is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of sugars present in the sample.

Linearity with total V / sample V = 50: up to 8 g/l

Precision: CV% < 1.5

**STEROGLASS S.r.l.**

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Calibration and control: standard SQPE063151 glucose-fructose-sucrose 20 g/l with sucrose 10 g/l, glucose 5g/l and fructose 5 g/l.

GLICER – GLYCERIN – SQPE060138 – 2 x 100 ml

Produced during glycerol-pyruvic fermentation, it greatly contributes to the characteristics of the wine giving harmonic taste sensation, body and roundness.

It would be appropriate to determine it after every re-fermentation.

Colorimetric method is suggested on wine, new UV is the best option for fruit juices and must.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GLIC R1 (GLICR)	2 x 100	large
GLIC BL	2 x100	optional

Method principle and characteristics

Glycerol takes part in a series of enzymatic reactions that produce H₂O₂ by means of glycerol kinase and glycerol phosphate oxidase. Hydrogen peroxide reacts with p-chlorophenol and 4-aminoantipyrine in the presence of peroxidase, forming a red-colored quinoneimine compound. The color intensity, measured at 520 nm, is proportional to the amount of glycerol present in the sample.

Linearity with total V / sample V =30: up to 12.5 g/l

Precision: CV% < 2

Calibration and control: standard SQPE053232 glycerol 10 g/l.

GLICUV – UV GLYCEROL – SQPE078691 – 50 ml

Produced during glycerol-pyruvic fermentation, it greatly contributes to the characteristics of the wine giving harmonic taste sensation, body and roundness. It would be appropriate to determine it after every re-fermentation.

Colorimetric method is suggested on wine, new UV is the best option for fruit juices and must.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
GLIC R1A (GUV1A)	1 x 20	small
GLIC R1B (GUV1B)	1 x 20	small
GLIC R2 (GUV2)	1 x10	small

Method principle and characteristics

In the presence of ATP, glycerol is converted into glycerol-3-phosphate by the enzyme glycerol kinase, with the resulting formation of ADP. The ADP that is formed reacts with phosphoenolpyruvate in the presence of pyruvate kinase to form ATP and pyruvate. In the presence of the enzyme LDH the pyruvate is reduced to lactate, oxidizing NADH. The decrease in absorbance measured at 340 nm, resulting from the oxidation of NADH, is proportional to the amount of glycerol present in the sample.

Linearity with total V / sample V =50: up to 6 g/l

Precision: CV% < 2

Calibration and control: standard SQPE053232 glycerol 10 g/l.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



GLUCON – GLUCONIC ACID – SQPE076314 – 1 x 50 ml

Organic acid from glucose by oxidation of the aldehyde function caused by an enzyme present on molds. The analysis is carried out on the raw material, on must and on wine. It is very important for importers and users of the grapes from different backgrounds. The enzymatic way is the analytical method faster and more precise.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
D-GLUC R1 (GLUC1)	1 x 40	large
D-GLUC R2 (GLUC2)	1 x 10	small
D-GLUC BL	1 x 50	do not use

Method principle and characteristics

In the presence of the enzyme gluconate kinase, D-gluconic acid is phosphorylated to D-gluconate-6-phosphate, with the consequential conversion of ATP into ADP. D-gluconate-6-phosphate, in the presence of NADP which is reduced to NADPH, is decarboxylated by the enzyme 6-phosphogluconate dehydrogenase and becomes ribulose 5-phosphate. The increase in absorbance due to the formation of NADPH, measured at 340 nm, is proportional to the amount of D-gluconic acid present in the sample.

Linearity with total V / sample V=50: up to 3 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE053234 gluconic acid 1 g/l.

K+ - POTASSIUM – SQPE056387 – 1 x 100 ml

Dominant element between cations present in wine.

The excess of potassium is removed through the tartaric stabilization and the determination is carried out before and after this process in order to verify the good performance.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
K R1 (POTAR)	1 x 100	large
K BL (POTAB)	1 x 100	large

Do not shake the reagents before use; if this happens, wait for them to settle.

Method principle and characteristics

The potassium ion reacts with sodium tetraphenylborate, forming a colloidal solution. Other components eliminate interference and stabilize the suspension. Consequently, it generates a certain turbidity that is proportional to the potassium concentration, which is read spectrophotometrically at 500 nm.

Linearity with total V / sample V=150: up to 1500 mg/l

Precision: CV% < 10

Calibration and control: standard SQPE056388 potassium 1500 mg/l.

Notes: as for all turbidimetric methods, it is recommended that multipoint calibrations be performed and the accuracy checked with standards at different concentrations. Given that the products of reactions are dirtying, the cuvette cleaning procedure should be carried out at the end of the analysis. Since this is the determination of the potassium ion, do not interfere by using potassium hydroxide (use NaOH) as a basic agent for cleaning.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

LLATAU – L-LACTIC ACID AUTO – SQPE074739 – 125 ml

L-lactic acid plays a fundamental role in the production of red wines. Significant quantities of lactic acid are found in wine (1.5 - 2.5 g/l) following the malolactic fermentation carried out by the lactic acid bacteria of the genus *Oenococcus*. The conversion of malic acid into lactic acid involves raising the pH, decreasing astringency and increasing the complexity of the aromas. Generally 0.64 g of lactic acid plus CO₂ are obtained from 1 g of malic acid. The starting of malolactic fermentation, and therefore the bacterial development, is strongly influenced by environmental factors, such as: pH, temperature, alcohol content, and content of energy and nitrogenous sources.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
L-LAT AUTO R1 (LATA1)	1 x 40	large
L-LAT AUTO R2 (LATA2)	1 x 10	small

Method principle and characteristics

L-lactate is converted into pyruvate by the enzyme L-lactate dehydrogenase in the presence of NAD. Given that the reaction equilibrium is clearly in favor of L-lactate, the pyruvate that is formed is progressively removed from the reaction environment by the activity of the enzyme D-glutamate pyruvate transaminase, which transforms it into D-alanine.

The NADH that is formed from the oxidation of D-lactate is in a stoichiometric concentration with the analyte L-lactic acid. Thus the absorbance is measured at a wavelength of 340 nm.

Linearity with total V / sample V = 100: up to 3 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 L-lactic acid 3 g/l; multilevel standard SQPE081638 L-lactic acid 3.0 g/l (level 1), 2.4 g/l (level 2), 1.8 g/l (level 3), 1.2 g/l (level 4), 0.6 g/l (level 5).

L-LATT – L-LACTIC ACID – SQPE059112 – 5 x 20 ml

L-lactic acid plays a fundamental role in the production of red wines. Significant quantities of lactic acid are found in wine (1.5 - 2.5 g/l) following the malolactic fermentation carried out by the lactic acid bacteria of the genus *Oenococcus*. The conversion of malic acid into lactic acid involves raising the pH, decreasing astringency and increasing the complexity of the aromas. Generally 0.64 g of lactic acid plus CO₂ are obtained from 1 g of malic acid. The starting of malolactic fermentation, and therefore the bacterial development, is strongly influenced by environmental factors, such as: pH, temperature, alcohol content, and content of energy and nitrogenous sources.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
L-LAT R1 (LLAT1)	5 x 20	small
L-LAT R2 (LLAT2)	1 x 2.5	cup
L-LAT BL	2 x 50	do not use

Method principle and characteristics

L-lactate is converted into pyruvate by the enzyme L-lactate dehydrogenase in the presence of NAD. Given that the reaction equilibrium is clearly in favor of L-lactate, the pyruvate that is formed is progressively removed from the reaction environment by the activity of the enzyme D-glutamate pyruvate transaminase, which transforms it into D-alanine.

The NADH that is formed from the oxidation of D-lactate is in a stoichiometric concentration with the analyte L-lactic acid. Thus the absorbance is measured at a wavelength of 340 nm.

Linearity with total V / sample V = 150: up to 3 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 L-lactic acid 3 g/l; multilevel standard SQPE081638 L-lactic acid 3.0 g/l (level 1), 2.4 g/l (level 2), 1.8 g/l (level 3), 1.2 g/l (level 4), 0.6 g/l (level 5).

MALAUT – L-MALIC ACID – SQPE068206 – 125 ml

Organic acid naturally present in fruit (in must and wine initial concentration of from about 1.5 g / l to 3.5). In wine production the initial monitoring can determine the knowledge of ripening course and subsequently, during and after the malo-lactic fermentation, to verify the end of this fermentation.

On apple juices it represents the most important organic acid.

20

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
L-MAL AUTO R1 (MALA1)	2 x 50	large
L-MAL AUTO R2 (MALA2)	1 x 25	small

Method principle and characteristics

L-malic acid is oxidized to oxaloacetate by the enzyme L-malate dehydrogenase in the presence of NAD, and therefore with the formation of NADH. Given that the reaction equilibrium is shifted towards L-malate, there is the presence of the enzyme glutamate-oxaloacetate transaminase, which subtracts the oxaloacetate from the reaction environment. The NADH that is formed is in a stoichiometric concentration with L-malic acid, and this increase in NADH is measured at 340 nm.

Linearity with total V / sample V = 100: up to 5 g/l

Precision: CV% < 1.5

Calibration and control: multiparameter standard SQPE053234 L-malic acid 5 g/l; multilevel standard SQPE081638 L-malic acid 5.0 g/l (level 1), 4.0 g/l (level 2), 3.0 g/l (level 3), 2.0 g/l (level 4), 1.0 g/l (level 5).

MALICO – L-MALIC ACID - SQPE053689 – 5 x 20 ml

Organic acid naturally present in fruit (in must and wine initial concentration of from about 1.5 g / l to 3.5). In wine production the initial monitoring can determine the knowledge of ripening course and subsequently, during and after the malo-lactic fermentation, to verify the end of this fermentation.

On apple juices it represents the most important organic acid.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
L-MAL R1 (MAL1)	5 x 20	small
L-MAL R2 (MAL2)	1 x 2.5	cup
L-MAL BL	2 x 50	do not use

Method principle and characteristics

L-malic acid is oxidized to oxaloacetate by the enzyme L-malate dehydrogenase in the presence of NAD, and therefore with the formation of NADH. Given that the reaction equilibrium is shifted towards L-malate, there is the presence of the enzyme glutamate-oxaloacetate transaminase, which subtracts the oxaloacetate from the reaction environment. The NADH that is formed is in a stoichiometric concentration with L-malic acid, and this increase in NADH is measured at 340 nm.

Linearity with total V / sample V = 100: up to 5 g/l

Precision: CV% < 1.5



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



Calibration and control: multiparameter standard SQPE053234 L-malic acid 5 g/l; multilevel standard SQPE081638 L-malic acid 5.0 g/l (level 1), 4.0 g/l (level 2), 3.0 g/l (level 3), 2.0 g/l (level 4), 1.0 g/l (level 5).

Mg²⁺ – MAGNESIUM – SQPE056389 – 2 x 100 ml

Magnesium is one of the cations most abundant in living cells, acting as a cofactor for over 300 enzymes and present in nature in different matrix, plays an important role in the precipitation of colloids. However, the concentration may be significantly lower due to certain amino acids and chelating agents that can complex magnesium, making it unavailable to microorganisms. Magnesium can be added to the must by adding magnesium sulfate.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
MAGN R1 (MGR)	1 x 100	small
MAGN R2 (MGR)	1 x 100	small
MAGN BL (MGB)	2 x 200	small

Mix reagents R1 and R2 in equal parts (1:1) to make up the MGR reagent. Remains stable for 2 days if stored at 4°C.

Method principle and characteristics

In an alkaline environment, calmagite binds with the magnesium present in the sample, forming a chromophore complex proportional to the quantity of analyte. Specific chelating agents eliminate interference due to the presence of other samples. Spectrophotometric measurement is done at 520 nm.

Linearity with total V / sample V = 100: up to 100 mg/l

Precision: CV% < 2

Calibration and control: standard SQPE056390 magnesium 24.3 mg/l.

NH₂ – AMINO NITROGEN - SQPE054974 – 2 X 60 ml

Sources of nitrogen are a fundamental factor for the multiplication and development of yeasts, the protagonists of fermentation. These sources are present in the must, but are often added, for example in the hottest years, to obtain a proper level of RAN (Readily Assimilable Nitrogen).

The α -amino nitrogen is identified in the amino acid fraction which is metabolized and catabolized by the yeasts at differing intensities according to the fermentation process. These microorganisms generally use the ammoniacal fraction in the early stages of fermentation, while in the later stages they mostly use the amino acid fraction.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
NH ₂ R1	2 X 50	large
NH ₂ R2	1 x 25	small
NH ₂ BL	2 x 60	do not use

Method principle and characteristics

In a basic environment with controlled pH, the amino group of free α -amino acids (NH₂) reacts with ortho-phthalaldehyde in the presence of N-acetylcysteine, forming chromophoric derivatives (isoindolic) in stoichiometric quantities with the free amino groups, which show an absorption band from 335 to 345 nm. The optical density read at 340 nm is proportional to the quantity of α -amino nitrogen present in the sample under examination.

Linearity with total V / sample V = 75: up to 200 mg/l

Precision: CV% < 2

Calibration and control: standard SQPE055107 α -amino nitrogen 160 mg/l



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

NH₃ – AMMONIACAL NITROGEN - SQPE054975 -2 X 50 ml

Sources of nitrogen are a fundamental factor for the multiplication and development of yeasts, the protagonists of fermentation. These sources are present in the must, but are often added, for example in the hottest years, to obtain a proper level of RAN (Readily Assimilable Nitrogen).

Ammoniacal nitrogen is normally metabolized by the yeast at differing intensities according to the fermentation.

These microorganisms generally use the ammoniacal fraction in the early stages of fermentation, while in the later stages they mostly use the amino acid fraction.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
NH3 R1	2 x 40	large
NH3 R2	1 x 20	small
NHH3 BL	2 x 50	do not use

Method principle and characteristics

Ammoniacal nitrogen reacts with α -ketoglutarate in the presence of NADH under the catalytic action of the enzyme glutamate dehydrogenase, producing NAD and L-glutamate.

The photometric measurement at 340 nm of the quantity of NADH oxidized to NAD indicates the quantity of ammoniacal nitrogen present in the sample.

Linearity with total V / sample V = 38: up to 100 mg/l

Precision: CV% < 2

Calibration and control: standard SQPE054993 ammoniacal nitrogen 50 mg/l.

PIRUV – PYRUVIC ACID – SQPE056391 – 5 x 20 ml

Organic acid produced by glycerol-pyruvic fermentation.

In fermentation on wine is not just a pure alcoholic fermentation, part of the sugar molecules are degraded by fermentation glycerol-pyruvic in glycerin and pyruvic acid. Their determination is performed during fermentation, especially because it combines highly SO₂.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
PIRU R1 (PIRR1)	5 x 20	small
PIRU R2 (PIRR2)	1 x 2.5	cup
PIRU BL	2 x 50	do not use

Method principle and characteristics

Pyruvic acid is reduced to D-lactic acid in the presence of the enzyme lactate dehydrogenase, with the consequential oxidation of NADH to NAD. The amount of NAD formed is proportional to the amount of pyruvic acid, and the decrease in absorbance, due to the disappearance of NADH, is measured at 340 nm.

Linearity with total V / sample V = 25: up to 0.6 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE056392 pyruvic acid 0.5 g/l.

Notes: to improve the performance of the analysis, it is recommended that structured red wines be decolorized with total polyphenols > 2500 mg/l.

POLIB – WHITE WINE POLYPHENOLS – SQPE054970 – 3 x 100 ml

Polyphenols accumulate mainly in the grape peel and seeds during the development of the berry and during the ripening of the fruit, therefore they are found in the must and especially in the wines undergoing maceration. The following classes of compounds fall within the category of polyphenols of oenological interest: hydrolysable tannins, condensed tannins, anthocyanins, flavonols and hydroxycinnamic acids. The concentration of total polyphenols depends on the grape variety, the environment, and technological factors.

These compounds determine a large part of the color and flavor quality of the wine, as they are at the basis of the astringency, the bitterness, the longevity and the evolutionary potential of the product. They are analyzed starting from the ripening of the grapes up to the bottling of the wine, given the numerous changes to which these compounds are subject.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
POLI3 R1 (POLI1)	2 x 100	large
POLI3 R2 (POLI2)	1 x 100	small

Method principle and characteristics

The potassium ion reacts with sodium tetraphenylborate, forming a colloidal solution. Other components eliminate interference and stabilize the suspension. Consequently, it generates a certain turbidity that is proportional to the potassium concentration, which is read spectrophotometrically at 500 nm.

Linearity with total V / sample V = 30: up to 600 mg/l

Precision: CV% < 10 in wine

Calibration and control: standard SQPE054995 polyphenols 3 g/l (anhydrous gallic acid).

Notes: This method is suitable for white wines with a polyphenol's concentration of up to 600 mg/l. To improve the accuracy of the method, it is advisable: to verify that the Blank value is less than 1; to perform the analysis of the polyphenols alone (not together with other methods); to clean the needle before the analysis.

POLIT – TOTAL POLYPHENOLS – SQPE054970 – 3 x 100 ml

They are extremely important because they characterize organoleptic aspect. Already present in fruit due to continually polymerization processes between anthocyanins and tannins evolve over time. It's important even to check polyphenols in fruit in order to determine the proper maturation for harvesting.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
POLI3 R1 (POLI1)	2 x 100	large
POLI3 R2 (POLI2)	1 x 100	small

Method principle and characteristics

In a strongly basic environment, the Folin-Ciocalteu reagent oxidizes the hydroxyl groups of the polyphenols present in the sample, so that the products of reduction formed are proportional to the analyte. These blue-colored products are detected spectrophotometrically at 700 nm.

Linearity with total V / sample V = 101: up to 3000 mg/l

Precision: CV% < 1.5

Calibration and control: standard SQPE054995 polyphenols 5 g/l.

Notes: This method is suitable for red wines with a polyphenol's concentration between 600 and 3000 mg/l. To improve the accuracy of the method, it is advisable: to verify that the Blank value is less than 1; to perform the analysis of the polyphenols alone (not together with other methods); to clean the needle before the analysis.

**PROT – PROTEIN – SQPE076312 – 1 x 50 ml**

The protein concentration in wine can vary from a few mg/l up to about 100 mg/l. These substances come from the grapes, and their content increases with the plant's response to stress in the case of disease. The oenological interest in proteins regards white wines, given the ability of proteins to impart turbidity to the product. In fact, the denaturation of proteins, such as that caused by heating, followed by their interaction with other compounds such as tannins, leads to the formation of invisible insoluble particles that aggregate and cause the wine to become cloudy. Generally, and when necessary, protein stability is achieved by means of physical treatments (ultrafiltration), the addition of tannins, proteolytic enzymes, polysaccharides, hyperoxygenation of the must or aging on the lees. Bentonite remains the most common additive for removing proteins.

To assess the degree of protein stability technicians, interpret the NTU values individually, whereas considering the protein concentration, the following stability thresholds can be considered, in reference to the type of product: below 25 mg/l, stable; 25 to 50 mg/l, at risk; and above 50 mg/l, unstable.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
PROT R1 (PROT1)	1 x 12.5	large
PROT R2 (PROT2)	1 x 50	small
PROT BL	1 x 25	do not use
STANDARD	1 x 10	Use as sample

As part of each analytical cycle, perform the analysis of two standards with a concentration of 50 mg/l (0.1 ml standard + 0.9 ml PROT BL) and 25 mg/l (1 ml standard at 50 mg/l + 1 ml of PROT BL), and then realign the results. Mix carefully to obtain a homogeneous solution, given the high viscosity of the standard.

The R1 reagent must be diluted 1:10 with distilled water (1 ml R1 + 9 ml H₂O), while the R2 reagent must be shaken gently before use.

Method principle and characteristics

The sample is given a pretreatment with a basic solution to favor the solubilization of the proteins, which then react with the Brilliant Blue G chromogen. The protein-chromogen complex generates a decrease in absorbance at 465 nm and an increase in absorbance at 595 nm. Therefore the 595 nm/465 nm ratio is a function of the concentration of proteins present in the sample.

Linearity with total V / sample V = 5.2: up to 100 mg/l

Precision: CV% < 5 in white wine

Calibration and control: standard proteins 1000 mg/l supplied with the Work Kit.

Notes: For samples containing more than 100 mg/l of protein, dilute the sample 1:4 with the BL reagent and multiply the result by 5. To increase the accuracy of the result, it is recommended that each analysis be performed three times and the result be considered as the average of the three values obtained.

**STEROGLASS S.r.l.**

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737



S – SUCROSE - SQPE063020 – 10 ml

In the F&B industry the sugar concentration is determined to establish the maturity of the fruit. In addition, the dosage of sugar for secondary fermentation is a common practice in the production. Sucrose may be present in traces, but it must be monitored in the production process

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
reagent R1 (GFAU1)	5 x 20	small
reagent R1B (GFAU2)	1 x 2.5	small
SACC	10	small

Reagents R1 and R1B are sold as GLUCOSE FRUCTOSE AUTO with the code SQPE068207.

Method principle and characteristics

Sucrose is hydrolyzed by the enzyme invertase into D-glucose and D-fructose. In the presence of phosphoglucose isomerase, D-fructose is converted into the D-glucose isomer, which is transformed into glucose-6-phosphate by the enzyme hexokinase, and is then oxidized to gluconate-6-phosphate by the enzyme glucose-6-phosphate dehydrogenase while reducing the NADP to NADPH. Thus the increase in absorbance, measured at 340 nm, due to the NADPH is proportional to the amount of glucose present in the sample.

Linearity with total V / sample V = 100: up to 6 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE063151 glucose-fructose-sucrose 20 g/l with sucrose 10 g/l, fructose 5g/l and glucose 5g/l.

Notes: This method is suitable for samples containing low concentrations of glucose and sucrose.

SO2L – FREE SULFUR DIOXIDE – SQPE056384 – 2 x 100 ml

Antioxidant and antiseptic functions of SO₂ are well known. SO₂ is added to several food compounds to preserve them. In wine the determination is carried out during the whole process of elaboration, from the must to the bottle. Free SO₂ is monitored especially during aging process.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
SO2L R1 (SO2LR)	2 x 80	small
SO2L R2 (SO2LB)	1 x 40	small
SO2L R3	1 x 5	0.25 parts in SO2LB

Both reagents must be prepared as follows: SO₂LR = 8 parts of R1 + 2 parts of R2 + 0.25 parts of distilled H₂O; SO₂LB = 8 parts of R1 + 2 parts of R2 + 0.25 parts of R3. The two reagents remain effective for 24 hours from the time of preparation.

Method principle and characteristics

The free sulfur dioxide in the sample reacts selectively with a chromogenic hydrochloride agent, developing a violet colored compound characterized by an absorption peak at 578 nm. Therefore the quantity of colored compound that is formed is proportional to the free SO₂ present in the sample.

Linearity with total V / sample V = 74: up to 100 mg/l

Precision: CV% < 3

Preparation of a standard of 80 mg/l of free sulfur dioxide: weigh 1.23g of sodium metabisulfite and put it in a 100ml volumetric flask. Add 40ml of distilled water, and shake until completely dissolved. Bring to volume with distilled water. Take 1 ml from this solution, put it in another 100 ml volumetric flask, and bring to volume with distilled water.



STEROGLASS S.r.l.

apparecchi in vetro soffiato
forniture per laboratori chimici

Cap. Soc. € 103.000,00 i.v.

Strada Romano di Sopra, 2/C - 06132 - S. Martino in Campo - Perugia

Tel. 075.609091 (r.a.) Fax 075.6090950

e-mail: info@sterglass.it ; sterglass.amm@pec.collabra.it ; internet: www.sterglass.it

R.I. / C.F. / P. IVA 01870870548 - VAT IT 01870870548 R.E.A. N. 164737

SO2T –TOTAL SULFUR DIOXIDE – SQPE060413 – 5 x 20 ml

Antioxidant and antiseptic functions of SO₂ are well known. SO₂ is added to several food compounds to preserve them. In wine the determination is carried out during the whole process of elaboration, from the must to the bottle. The official method of distillation is not practical; the usual one, Ripper, gives good results on whites but not on red wines due to the interference of polyphenols, tannins and color.

The total SO₂ is the sum of the different fractions present and is considered in reference to the legal limits. If a high total SO₂ value is required to obtain effective free values, this indicates that a large part of the sulfur dioxide has combined with specific compounds

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
SO2T R1 (SO2T1)	2 x 80	small
SO2T R2 (SO2T2)	1 x 40	small
SO2T BL	1 x 5	do not use

Method principle and characteristics

The total sulfur dioxide in the sample reacts selectively with a chromogenic agent having a disulfide bridge, giving rise to a yellow colored compound characterized by an absorption peak at 420 nm. The intensity of the color formed is proportional to the amount of total sulfur dioxide present in the sample.

Linearity with total V / sample V = 50: up to 300 mg/l

Precision: CV% < 2

Preparation of a standard of 160 mg/l of free sulfur dioxide: weigh 1.23g of sodium metabisulfite and put it in a 100ml volumetric flask. Add 40ml of distilled water, and shake until completely dissolved. Bring to volume with distilled water. Take 1 ml from this solution, put it in a 50 ml volumetric flask, and bring to volume with distilled water.

TARTAR – TARTARIC ACID – SQPE070208 – 179 ml

It is the organic acid is specific of grapes, the most important of fixed acids, the stronger, the more dissociated, and the most resistant to the action of the decomposing bacteria. In wine its concentration decreases by precipitation of potassium bitartrate during the process of tartaric stabilization; it is determined on the must, on the wine after tartaric stabilization and before bottling.

COMPONENTS SUPPLIED (ID abbrev.)	V(ml)	BOTTLE-REAGENT
TART R1 (TART1)	1 x 136	large
TART R2 (TART2)	1 x 34	small
TART BL	2 x 85	do not use

Method principle and characteristics

Tartaric acid reacts with vanadic acid in an acid medium to develop a colored complex, the concentration of which is directly proportional to the amount of tartaric acid in the sample. The absorbance is measured at 492 nm.

Linearity with total V / sample V = 22: up to 6 g/l

Precision: CV% < 1.5

Calibration and control: standard SQPE054994 tartaric acid 5 g/l.

Notes: to improve the performance of the analysis, it is recommended that structured red wines be decolorized with total polyphenols > 2500 mg/l. In addition, the result found should be added to +0.14 (HM).