

Redox Electrode for pH 7 and pH 70

5MT cable+ BNC

Item Code: 32105102



DATA SHEET

Features

Features	Maintenance free
Connection	BNC
	Combined Redox electrode. General purpose plastic body. 1 m fixed cable, BNC plug.
	Gel filling with porous septum. Maintenance-free.
	Range $\pm 1,000$ mV, Temperature 0...60 °C. Dimensions (WxD) 120 x 12 mm
Description	
Diaphragm	porous septum
Body diameter	12 mm
Electrolyte	Gel
Image 2	https://app.onpage.it/api/storage/a56538eba66645f81fa28bc724f3e5b83e437
Body lenght	120 mm
Cable length	5 m
Body material	Plastic
Name	REDOX / 5
mV ORP measuring range	-1000 ... +1000 mV