

pH Electrode 201T DHS for pH 7 / pH 70 / pH50

Item Code: 32200103



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
Description	Digital combined pH electrode with DHS technology. General purpose plastic body. Gel filling with porous septum. Maintenance-free. Dimensions (WxD) mm 120x12. pH 0...14 - Temperature 0...60 °C
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	Gel
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Plastic
Name	201 T DHS
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Ag / AgCl