

# XS Sensor 2-Pore T DHS NTC with 1mt cable BNC

Item Code: 32200113



## DATA SHEET

### Features

Application	Ideal for meat and cheese
Cable	Si
Connection	Fixed cable / BNC
Description	Digital combined pH electrode with DHS technology. Glass body. Maintenance-free polymer filling. The 6 mm tip makes it ideal for meat and cheese.
	Dimensions (WxD) mm 35x6 - pH 0...14, temperature 0...60 °C
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	<a href="https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6">https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6</a>
Body lenght	95 / 35 mm
Cable length	1 m
Body material	Glass
Name	2 PORE DHS
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Covered with Ag