

XS Sensor Standard T DHS BNC

Item Code: 32200123



DATA SHEET

Features

Cable	Si
Certificate	Yes
Connection	Fixed cable / BNC
Digital combined pH electrode with DHS technology. Cylindrical membrane, general purpose. Refillable 3M KCl electrolyte. Excellent chemical resistance and easy to clean. With serial number and certificate.	
Dimensions (WxD) mm 120x12 - pH 0...14 Temperature - 10...100 °C	
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	STANDARD DHS
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Ag / AgCl Encapsulated