

# XS Sensor 2 Pore F Temp-BNC DHS

Item Code: 32200143



## DATA SHEET

### Features

Application	Ideal for meat and cheese and penetration measures
Cable	Si
Connection	Fixed cable / BNC
Digital combined pH electrode with DHS technology. Plastic body(PEEK) and glass membrane. Polymer-filled, maintenance-free. Thanks to the 6 mm tip it is ideal for meat and cheese and penetration measurements.	
Dimensions (WxD) mm 35x6 - pH 0...14, temperature 0...60 °C	
Description	
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	<a href="https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6">https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6</a>
Body lenght	90 / 50 mm
Cable length	1 m
Body material	POM / Glass
Name	2 PORE F DHS
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Covered with Ag