

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Gel S7

Item Code: 32200203



DATA SHEET

Features

Connection	Screw head S7
	Combined pH electrode, screw head S7. General purpose glass body. Cylindrical membrane. Maintenance-free gel filling
	Dimensions (WxD) mm 120x12 - pH 0...14 -Temperature -10...60 °C
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	Gel
Image 2	https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC
Body lenght	120 mm
Body material	Glass
Name	Gel
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +60 °C
Reference system	Ag /AgCl