

XS Sensor Micro S7

Item Code: 32200213



DATA SHEET

Features

Application	For small sample volumes < 2 ml
Connection	Screw head S7
	Combined pH electrode, screw head S7. Micro 3" glass body. 3M KCl refillable electrolyte. For small sample volumes <2 ml.
	Dimensions (WxD) mm 60 x 3 - pH 0...14 - Temperature 0...100 °C
Description	
Diaphragm	Ceramic
Body diameter	12 / 3 tip mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/61e6117613356a6967fd79aea59e63c8e61e6c
Body lenght	120 / 60 mm
Body material	Glass
Name	Micro
pH measuring range	Vetro pH
temperature measurement range	0 ... +100 °C
Reference system	Ag / AgCl Encapsulated