

# XS Sensor Polymer S7

Item Code: 32200223



## DATA SHEET

### Features

Connection	Screw head S7
	Combined pH electrode, screw head S7. General purpose glass body. Cylindrical membrane. Polymer filling without porous septum. For fast response with high accuracy in samples that may contaminate normal diaphragms. Maintenance-free.
Description	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature - 10...80 °C
Diaphragm	Open hole
Body diameter	12 mm
Electrolyte	Polymer
Image 2	<a href="https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC">https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC</a>
Body lenght	120 mm
Body material	Glass
Name	Polymer
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
Reference system	Ag / AgCl Encapsulated