

XS Sensor Standard ORP S7

Item Code: 32200233



DATA SHEET

Features

Application	Redox measurements
Connection	Screw head S7
	Combined Redox electrode, screw head S7. Glass body. Filled with 3M KCL
Description	Dimensions (WxD) mm 120x12 - Temperature -10...100 °C
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Body lenght	120 mm
Body material	Glass
Name	Standard RFO
Number of diaphragm	3
mV ORP measuring range	-2000 ... +2000 mV
temperature measurement range	-10 ... +100 °C
Temperature sensor	-10 ... +100
Reference system	Ag / AgCl Encapsulated