

XS Sensor Food Slim S7

Item Code: 32200243



DATA SHEET

Features

Application	Penetration measurements for protein-containing samples
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body. Electrolyte Glycerin KCl. Suitable for protein-containing samples. 6 mm tip for penetration measurements.
	Dimensions (WxD) mm 25x6 - pH 0...14 - Temperature - 10...100 °C
Description	
Diaphragm	Ceramic
Body diameter	12 / 6 tip mm
Electrolyte	Glycerin KCl
Image 2	https://app.onpage.it/api/storage/801ca24af9c53185beb9809f247b26746903c
Body lenght	95 /25 mm
Body material	Glass
Name	Food Slim
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated