

# XS Sensor Micro P S7

Item Code: 32200253



## DATA SHEET

### Features

|                               |                                                                                                                                                                                              |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application                   | Small volumes in biochemical and microbiological applications                                                                                                                                |
| Connection                    | Screw head S7                                                                                                                                                                                |
|                               | Combined pH electrode, screw head S7. Micro 3" glass body. With refillable Glycerin KCl electrolyte. For measurements in small sample volumes, biochemical and microbiological applications. |
|                               | Dimensions (WxD) mm 60 x 3 - pH 0...14 - Temperature 0...100 °C                                                                                                                              |
| Description                   |                                                                                                                                                                                              |
| Diaphragm                     | Ceramic                                                                                                                                                                                      |
| Body diameter                 | 12 / 3 tip mm                                                                                                                                                                                |
| Electrolyte                   | Glycerin KCl                                                                                                                                                                                 |
| Image 2                       | <a href="https://app.onpage.it/api/storage/b5ad038a4fc5a7fa270ddff64d592fb6067bcf">https://app.onpage.it/api/storage/b5ad038a4fc5a7fa270ddff64d592fb6067bcf</a>                              |
| Body lenght                   | 120 / 60 mm                                                                                                                                                                                  |
| Body material                 | Glass                                                                                                                                                                                        |
| Name                          | Micro P                                                                                                                                                                                      |
| pH measuring range            | Vetro pH                                                                                                                                                                                     |
| temperature measurement range | 0 ... +100 °C                                                                                                                                                                                |
| Reference system              | Ag / AgCl Encapsulated                                                                                                                                                                       |