

XS Sensor Standard + Cable BNC

Item Code: 32200263



DATA SHEET

Features

Cable	Si
Certificate	Yes
Connection	Fixed cable / BNC
Combined pH electrode, screw head S7. Glass body, general purpose. Refillable 3M KCl electrolyte. Excellent chemical resistance and easy to clean. With serial number and certificate.	
Dimensions (WxD mm) 120x 12 - pH 0...14 Temperature - 10...100 °C	
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/b1b66d3034c2f6843aa534656db82fcd645
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	Standard
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated