

XS Sensor 2-Pore Slim S7

Item Code: 32200283



DATA SHEET

Features

Application	Ideal for meat and cheese and penetration measures
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body with double open hole. Polymer filling, maintenance-free. 4.5 mm tip: ideal for meat and cheese and penetration measurements.
	Dimensions (WxD) mm 35x4,5 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Open hole
Body diameter	12 / 4.5 tip mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/801ca24af9c53185beb9809f247b26746903c
Body lenght	100 / 35 mm
Body material	Glass
Name	2 Pore Slim
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Encapsulated ag