

XS Sensor 2-Pore K with Stainless steel knife S7

Item Code: 32200303



DATA SHEET

Features

Application	Ideal for meat and penetration measurements
Connection	Screw head S7
	Combined pH electrode, screw head S7. Steel body with interchangeable blade. Polymer filling, maintenance-free.
	The 6 mm tip makes it ideal for meat. For penetration measurements.
	Dimensions (WxD) mm 35x6 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/b9fc9dbf5f3fd2e2e8c8c1ffb72e0ff57be3db1
Body lenght	110 / 35 mm
Body material	Glass / Steel
Name	2 Pore K
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Encapsulated ag