

XS Sensor 2-Pore NTC with 1mt cable BNC

Item Code: 32200323



DATA SHEET

Features

Application	Ideal for meat and cheese
Cable	Si
Connection	BNC + RCA
Combined pH electrode, screw head S7. Glass body with double open hole. Polymer filling, maintenance-free. 6 mm tip: ideal for meat and cheese and penetration measurements.	
Dimensions (WxD)mm 35x6 - pH 0...14 - Temperature 0...60 °C	
Description	
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	95 / 35 mm
Cable length	1 m
Body material	Glass
Name	2 Pore T
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Encapsulated ag