

Via della Meccanica, 25, 41012 Carpi (MO) - Italy  
+39 059 653274 - info@giorgiobormac.com

# XS Sensor 2-Pore S7

Item Code: 32200333



## DATA SHEET

### Features

Application	Ideal for meat and cheese and penetration measures
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body with double open hole. Polymer filling, maintenance-free. 6 mm tip: ideal for meat and cheese and penetration measurements.
Description	Dimensions (WxD) mm 35x6 - pH 0...14 - Temperature 0...60 °C
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	<a href="https://app.onpage.it/api/storage/801ca24af9c53185beb9809f247b26746903c">https://app.onpage.it/api/storage/801ca24af9c53185beb9809f247b26746903c</a>
Body lenght	100 / 35 mm
Body material	Glass
Name	2 Pore
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Encapsulated ag