

XS Sensor 2-Pore K without knife

Item Code: 32200343



DATA SHEET

Features

Application	Ideal for meat and penetration measurements
Connection	Screw head S7
	Combined pH electrode, S7 screw head. Steel body, with thread to accommodate the blade. Polymer filling, maintenance-free. The 6 mm tip is ideal for meat. For penetration measurements.
	Dimensions (WxD) mm 35x6 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/b9fc9dbf5f3fd2e2e8c8c1ffb72e0ff57be3db1
Body lenght	110 / 35 mm
Body material	Glass / Steel
Name	2 Pore K
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Encapsulated ag