

Via della Meccanica, 25, 41012 Carpi (MO) - Italy  
+39 059 653274 - info@giorgiobormac.com

# Knife for 2-Pore

Item Code: 32200353



## DATA SHEET

### Features

Combined pH electrode, S7 screw head. Steel body, with thread to accommodate the blade.

Polymer filling, maintenance-free.

The 6 mm tip is ideal for meat. For penetration measurements.

Dimensions (WxD)mm 35x6 - pH 0...14 - Temperature 0...60 °C

Description

Image 2

<https://app.onpage.it/api/storage/4b1d019056e5a336d7ded7f5a5bf653726e59>

Name

Blade