

XS Sensor Standard S7

Item Code: 32200363



DATA SHEET

Features

Certificate	Yes
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body, general purpose. Refillable 3M KCl electrolyte. Excellent chemical resistance and easy to clean. With serial number and certificate
Description	Dimensions (WxD mm) 120x 12 - pH 0...14 Temperature - 10...100 °C
Diaphragm	Ceramics
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC
Body lenght	120 mm
Body material	Glass
Name	Standard
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated