

# XS Sensor Flow S7

Item Code: 32200373



## DATA SHEET

### Features

Application	Suspension measurements and low conductivity samples
Certificate	Yes
Connection	Screw head S7
Combined pH electrode, screw head S7. Glass body. Refillable 3M KCl electrolyte. Collar diaphragm for measurements in suspensions and low conductivity samples. With serial number and certificate.	
Description	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature - 10...80 °C
Diaphragm	Teflon collar
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	<a href="https://app.onpage.it/api/storage/44bd457369c9068464ac2d0a697a161ee750t">https://app.onpage.it/api/storage/44bd457369c9068464ac2d0a697a161ee750t</a>
Body lenght	120 mm
Body material	Glass
Name	Flow
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
Reference system	Ag / AgCl Encapsulated