

XS Sensor Semi Micro S7

Item Code: 32200383



DATA SHEET

Features

Application	Measurements in test tubes
Connection	Screw head S7
	Combined pH electrode, screw head S7. Micro 6" glass body. KCl 3M electrolyte. Suitable for measurements in test tubes.
	Dimensions (WxD) mm 100 x 6 - pH 0...14 - Temperature 0...100 °C
Description	
Diaphragm	Ceramic
Body diameter	12 / 6 tip mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/d09c2a4b66b100afe43c3868b3515703eaffd7
Body lenght	150 / 100 mm
Body material	Glass
Name	Micro Seeds
pH measuring range	Vetro pH
temperature measurement range	0 ... +100 °C
Reference system	Ag / AgCl Encapsulated