

XS Sensor Zava trode S7

Item Code: 32200403



DATA SHEET

Features

Application	Penetration measurements for protein-containing samples
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body. Electrolyte Glycerin KCl. Suitable for protein-containing samples. 4 mm tip for penetration measurements.
	Dimensions (WxD)mm 120x4 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Ceramic
Body diameter	12 / 4 tip mm
Electrolyte	Glycerin KCl
Image 2	https://app.onpage.it/api/storage/ba6915901be2ae99391363ad734c1c9580907
Body lenght	190 / 120 mm
Body material	Glass
Name	Zavatrode
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Ag /AgCl