

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Flat

Item Code: 32200413



DATA SHEET

Features

Application	Surface measurements
Connection	Screw head S7
	Combined pH electrode, screw head S7. Glass body. Polymer filling. With flat membrane for surface measurement.
	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Teflon
Body diameter	12 mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/46a41b2e36ca7209394b0333a6ecb246a364
Body lenght	120 mm
Body material	Glass
Name	Flat
Number of diaphragm	3
pH measuring range	0 ... 14 pH
temperature measurement range	0 ...+ 60 °C
Reference system	Encapsulated ag