

XS Sensor Standard HA S7

Item Code: 32200423



DATA SHEET

Features

Application	Highly alkaline products
Certificate	Yes
Connection	Screw head S7
Low error combined pH electrode in alkaline environment, screw head S7. Glass body and cylindrical membrane. Refillable 3M KCl electrolyte. Excellent chemical resistance and easy to clean. With serial number and certificate.	
Dimensions (WxD mm) 120x 12 - pH 0...14 Temperature - 10...100 °C	
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC
Body lenght	120 mm
Body material	Glass
Name	Standard HA
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated