

# XS Sensor Standard HF S7

Item Code: 32200433



## DATA SHEET

### Features

Application	Products with HF
Certificate	Yes
Connection	Screw head S7
Combined pH electrode, screw head S7. Glass body. With special membrane for samples containing hydrofluoric acid (0.01M/200 mg at 20°C). Refillable 3M KCl electrolyte. Excellent chemical resistance and easy to clean.	
Dimensions (WxD mm) 120x 12 - pH 0...14 Temperature - 10...100 °C	
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	<a href="https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC">https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC</a>
Body lenght	120 mm
Body material	Glass
Name	Standard HF
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated