

XS Sensor MICRO-Special S7

Item Code: 32200453



DATA SHEET

Features

Application	NMR tube measurements
Connection	Screw head S7
	Combined pH electrode, screw head S7. Micro 3" glass body. 3M KCl electrolyte. For sample volumes below 2 ml. Ideal for measurements in NMR tubes.
	Dimensions (WxD) mm 180x3 - pH 0...14 - Temperature 0...100 °C
Description	
Diaphragm	Ceramic
Body diameter	12 / 3 tip mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/411913d7a34fa95ea05d0eebc6f3cf65af126d
Body lenght	220 / 180 mm
Body material	Glass
Name	Micro Special
pH measuring range	Vetro pH
temperature measurement range	0 ... +100 °C
Reference system	Ag / AgCl Encapsulated