

XS Sensor Standard T BNC

Item Code: 32200473



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
	Combination pH electrode, S7 screw head. Glass body, general purpose. 3M refillable KCl electrolyte. Excellent chemical resistance and easy to clean.
	Supplied without certificate.
Description	Dimensions (WxØ mm) 120x 12 - pH 0...14 Temperature - 10...100 °C
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	Standard T
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Temperature sensor	Si
Reference system	Ag / AgCl Encapsulated