

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Standard Ag S7

Item Code: 32200483



DATA SHEET

Features

Application	Argentometric titrations
Connection	Screw head S7
	S7 screw head. Glass body. KNO3 1M filler.
	Electrode for argentometric titrations.
Description	Dimensions (WxD) mm 120X12- Temperature 0...80 °C
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KNO3
Body lenght	120 mm
Body material	Glass
Name	Standard AG
mV ORP measuring range	-2000 ... +2000 mV
temperature measurement range	-10 ... +80 °C
Temperature sensor	-10 ... +80
Reference system	Ag / AgCl Encapsulated