

XS Sensor Polymer PLAST

Item Code: 32200493



DATA SHEET

Features

Combined pH electrode, S7 screw head, with unbreakable epoxy body. Polymer filling one open porous septum and one ceramic septum for fast response with high accuracy in samples that may contaminate normal diaphragms.
Maintenance-free

Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C

| | |
|-------------------------------|---|
| Description | |
| Diaphragm | Ceramic + Open Hole |
| Body diameter | 12 mm |
| Electrolyte | Polymer |
| Image 2 | https://app.onpage.it/api/storage/932ba8eb63be36a4e024184d4cf4bdd747777 |
| Body lenght | 120 mm |
| Body material | Plastic |
| Name | Polymer Plast |
| Number of diaphragm | 2 |
| pH measuring range | 0 ... 14 pH |
| temperature measurement range | 0 ... +60 °C |
| Reference system | Ag / AgCl Encapsulated |