

XS Sensor Polymer PLAST BNC

Item Code: 32200503



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
Description	Combined pH electrode, S7 screw head, with unbreakable epoxy body. Polymer filling one open porous septum and one ceramic septum for fast response with high accuracy in samples that may contaminate normal diaphragms. Maintenance-free. Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C
Diaphragm	Ceramic + Open Hole
Body diameter	12 mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/b1b66d3034c2f6843aa534656db82fdbcdb649
Body lenght	120 mm
Cable length	1 m
Body material	Plastic
Name	Polymer Plast
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Ag / AgCl Encapsulated