

XS Sensor Polymer PLAST DIN

Item Code: 32200513



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / DIN
	Combined pH electrode, S7 screw head, with unbreakable epoxy body. Polymer filling one open porous septum and one ceramic septum for fast response with high accuracy in samples that may contaminate normal diaphragms.
	Maintenance-free
	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Ceramic + Open Hole
Body diameter	12 mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/a4b2557c4a3540c981625b1398761e6e5e319
Body lenght	120 mm
Cable length	1 m
Body material	Plastic
Name	Polymer Plast
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Reference system	Ag / AgCl Encapsulated