

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Flow Ag

Item Code: 32200543



DATA SHEET

Features

Application	Argentometric applications
Connection	Screw head S7
	S7 screw head. Glass body. Teflon diaphragm
	Electrode for argentometric titrations.
Description	Dimensions (WxD) mm 120X12- Temperature 0...80 °C
Diaphragm	Collar
Body diameter	12 mm
Body lenght	120 mm
Body material	Glass
Name	Flow AG
mV ORP measuring range	-2000 ... +2000 mV
temperature measurement range	0 ... +80 °C
Temperature sensor	0 ... +80