

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor 2-Pore steel S7

Item Code: 32200653



DATA SHEET

Features

Application	Ideal for meat and penetration measurements
Connection	Screw head S7
	Steel body. S7 head. Polymer filling, maintenance-free. 6 mm tip for penetration measurements. Ideal for meat.
	Dimensions (WxD) mm 50x6 - pH 0...14 - Temperature 0...80 °C
Description	
Diaphragm	Open hole
Body diameter	12 / 6 tip mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/5728df63ff1acee13f864f53440e02f244060c
Body lenght	105 / 50 mm
Body material	Glass / Steel
Name	2 Pore Steel
Number of diaphragm	2
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... 80 °C
Reference system	Encapsulated ag