

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

Sensor Poly Plast ORP S7

Item Code: 32200663



DATA SHEET

Features

	S7 screw head, epoxy body.
	Maintenance-free gel filling. Porous ceramic septum.
Description	Dimensions (WxD) mm 120X12 - Temperature 0...60 °C
Image 2	https://app.onpage.it/api/storage/e6cbc36ea235d7ca3e8e94048fa7336c2ada
Name	Plast ORP