

XS Sensor Flow T BNC NTC 30K

Item Code: 32200723



DATA SHEET

Features

Application	Suspension measurements and low conductivity samples
Cable	Si
Certificate	Yes
Connection	Fixed cable / BNC
Description	Combined pH electrode, screw head S7. Glass body. Refillable 3M KCl electrolyte. Collar diaphragm for measurements in suspensions and low conductivity samples. With serial number and certificate.
	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature - 10...80 °C
	Teflon collar
Diaphragm	Teflon collar
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	Flow T
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Ag / AgCl Encapsulated