

XS Sensor Flow T DIN PT 1000

Item Code: 32200733



DATA SHEET

Features

Application	Suspension measurements and low conductivity samples
Cable	Si
Certificate	Yes
Connection	Fixed cable / DIN
Description	Combined pH electrode, screw head S7. Glass body. Refillable 3M KCl electrolyte. Collar diaphragm for measurements in suspensions and low conductivity samples. With serial number and certificate. Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature - 10...80 °C
Diaphragm	Teflon collar
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/15821a121b58514827a0d552259f10accd505f
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	Flow T
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
Temperature sensor	Si / spina banana 4 mm
Reference system	Ag / AgCl Encapsulated