

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor Food T

Item Code: 32200753



DATA SHEET

Features

Application	Dairy and food sector
Connection	BNC + RCA
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	Glycerin KCL
Image 2	Food T
Body lenght	120 mm
Body material	Glass
Name	Food T
Number of diapraghm	3
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... + 100 °C
Reference system	Ag / AgCl Encapsulated