

XS Sensor Electrode G 210 T Pt 1000 (BNC+Banana 4 mm) Plastic Body

Item Code: 32200803



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
	Combined pH electrode with built-in NTC temperature sensor. General purpose plastic body. Gel filling with porous septum. Maintenance-free. 1 m fixed cable with BNC+RCA plug
	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	Gel
Image 2	https://app.onpage.it/api/storage/125e9943f52137392210bbe758328ccd15e6e!
Body lenght	120 mm
Cable length	1 m
Body material	Plastic
Name	G 210
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina banana 4 mm
Reference system	Ag / AgCl