

XS Sensor Electrode G 214 T Pt 1000 (BNC+Banana 4 mm) Glass Body

Item Code: 32200813



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
	Combination pH electrode, S7 screw head. Glass body, general purpose. 3M refillable KCl electrolyte. Excellent chemical resistance and easy to clean.
	Supplied without certificate.
	Dimensions (WxØ mm) 120x 12 - pH 0...14 Temperature - 10...100 °C
Description	
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	https://app.onpage.it/api/storage/125e9943f52137392210bbe758328ccd15e6e!
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	G 214
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Temperature sensor	Si / spina banana 4 mm
Reference system	Ag / AgCl Encapsulated