

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

XS Sensor G 221 T

Item Code: 32200863



DATA SHEET

Features

Cable	Si
Certificate	Yes
Connection	Fixed cable / BNC
Diaphragm	Teflon collar
Body diameter	12 mm
Electrolyte	KCl 3M
Image 2	G 221 T
Cable length	1 m
Body material	Glass
Name	G 221 T
Serial number	Si
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +100 °C
Reference system	Ag / AgCl Encapsulated