

MULTISPIN 12 - Mini Microcentrifuge

6x2ml, 12.000rpm

Item Code: 40005012



DATA SHEET

Features

Power supply	220, 50-60 Hz V
Capacity max. (fixed angle rotor)	6 x 2 ml
Consumption	30 W
Dimension	176 x 196 x 122 mm
Warranty	2 anni
RPM settings	Yes
Name	MULTISPIN 12
Weight	1.5 kg
RCF max (fixed angle rotor)	7900 x g
Noise	45 dB
Control system	Microprocessor and LED display
Timer	1 - 99 min hh:mm
Type of motor	Brushless induction
maximum speed	12000 rpm