

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

WATERBATH ARGOLAB WB 22 Lt (with pump)

Item Code: 41101612



DATA SHEET

Features

Power supply	230 V
Min. height (lid closed)	150 mm
Capacity	22 L
Security class	2
Dimension	680 x 365 x 390 mm
Internal dimensions	500 x 295 x 200 mm
Dimension (plate)	450 x 265 mm
Display	LCD
Frequency	50 - 60 Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/b010929ff428c8ef6d25c0e9b63342f800280
Working method	forced circulation
Name	WB 22
Regulations	2006/95/EC - 2004/108/EC - EN61326-1:2006 - EN61010-1:2010 - EN61010-2-010:2003
Temperature homogeneity at 37 °C	±0,2 °C
Weight	19 kg
Recirculation pump	Yes
Pump flow rate max.	12 l/min
Absorbed power	1100 W
Tub width	200 mm
Resolution	0,1
Temperature resolution	0,1 °C
Overtemperature security	Electronic
Working temperature	+5 Ta...85 °C
Timer	99 : 59 hh:mm
Temperature variation at 37 °C	±0,1 °C
Useful volume	22 L