


Shaker Incubator SKI 8

Item Code: 41102022



DATA SHEET

Features

Accuracy (stirrer speed)	1 rpm
Power supply	230 V
Capacity	8 L
Heating accuracy control	0,1 °C
Plate dimension	450 x 450 mm
Dimension	610 x 840 x 540 mm
Stir excursion	20 mm
Frequency	50 / 60 Hz Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/4fbbfb63204ff27b9f27acbe1f3e9619ec7a1fd
Temperature indication	Yes
Time indication	Yes
Speed indication	Yes
Plate material	Metal
Stir movement	Orbital
Name	SKI 8
Regulations	2006/95/EC -2004/108/EC - EN61326-1:2006 - EN61010-1:2010 - EN61010-2-010:2003
Weight	55 kg
Absorbed power	1000 W
Heating power	1000 W
Over Temperature Protection	Yes
Heating Temperature Range	+ 5 Ta ... 60 °C
Timer	99:59 o continuo hh:mm
Stirring speed	40 ... 300 rpm