

Digital ultrasonic bath Mod. DU-100

Item Code: 41300363



DATA SHEET

Features

Capacity	10 L
Features	Bright LCD display showing current temperature, set temperature, timer, ultrasonic power level and operating frequency
Power control (%)	Yes
Tub dimension	300 x 240 x 150 mm
External dimensions	325 x 265 x 310 mm
Supplied as standard	Stainless steel basket and lid
Frequency	40.000/28.000 Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/138018a8dee43a12eb611b6e03a9a7f32f7e3e
Name	DU-100S
Regulations	2014/30/UE ; 2014/35/UE ; EN 61010-1: 2010 ; EN 61326-1: 2013 ; EN 61000-3-3: 2013 ; EN 61000-3-2: 2014
Number of transducer	4
Weight	8,2 kg
Heating power	500 W
Ultrasonic power	240 W
Set temperature	+5°C Ta ... +60 °C
Timer	00:99 hh:mm
Discharge valve	Yes