

Analogue ultrasonic bath Mod. AU-450

Item Code: 41300423



DATA SHEET

Features

	Analogue ultrasonic bath
	Temperature control up to 80°C
	Timer can be set from 1 to 30 minutes and infinite
Other features	ultrasonic frequency 40 kHz
Capacity	45 L
Power control (%)	No
Tub dimension	500 x 300 x 300 mm
External dimensions	530 x 325 x 500 mm
Supplied as standard	Stainless steel lid and basket
Frequency	40 Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/6aca51874163652b590758e28a329b7db1395
Name	AU-450
Regulations	2014/30/UE ; 2014/35/UE ; EN 61010-1: 2010 ; EN 61326-1: 2013 ; EN 61000-3-3: 2013 ; EN 61000-3-2: 2014
Number of transducer	12
Weight	17 kg
Heating power	1000 W
Ultrasonic power	750 W
Set temperature	+5°C Ta ... +80°C °C
Timer	00:30 e in continuo hh:mm
Discharge valve	Yes