

Analogue ultrasonic bath Mod. AU-220

Item Code: 41300433



DATA SHEET

Features

	Analogue ultrasonic bath
	Temperature control up to 80°C
	Timer can be set from 1 to 30 minutes and infinite
Other features	ultrasonic frequency 40 kHz
Capacity	22 L
Power control (%)	No
Tub dimension	500 x 300 x 150 mm
External dimensions	500 x 330 x 340 mm
Supplied as standard	Stainless steel lid and basket
Frequency	40 Hz
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/6aca51874163652b590758e28a329b7db1395
Name	AU-220
Regulations	2014/30/UE ; 2014/35/UE ; EN 61010-1: 2010 ; EN 61326-1: 2013 ; EN 61000-3-3: 2013 ; EN 61000-3-2: 2014
Number of transducer	8
Weight	13,2 kg
Heating power	500 W
Ultrasonic power	480 W
Set temperature	+5°C Ta ... +80°C °C
Timer	00:30 e in continuo hh:mm
Discharge valve	Yes