

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

ArgoLab Refrigerator FR 110 - blind door

Item Code: 41500012



DATA SHEET

Features

Power supply	220 - 240 V
Maximum shelf load	15 kg
Consumption	584 KW/year
Energy consumption (kWh/24h)	1.5 kWh/24h
Through-hole diameter	10 mm
External dimensions	600 x 560 x 805 mm
Internal dimensions	500 x 420 x 570 mm
Display	microprocessor/display
Minimum distance between shelves	13 mm
Defrosting function	Automatic and manual
Gradient at 5°C	1.3 °C
Image 2	https://app.onpage.it/api/storage/327a41c257bd4d0d7f14768c04c4d3ca1d47a1
Internal light	Yes
Name	FR 110
Number of shelves	3
Weight	50 kg
Power	120 W
temperature measurement range	2 ... 8 °C
Refrigeration	Dynamics
Noise	55 dB
Key Lock	Yes
Recovery time 5°C	12 m
Door type	Blind
Temperature variation at 5 °C	1,3 °C
Useful volume	110 L