

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

pH Electrode Model 201T + ATC 5MT

Item Code: 50002022



DATA SHEET

Features

Cable	Si
Connection	Fixed cable / BNC
	Combined pH electrode with built-in NTC temperature sensor. General purpose plastic body. Gel filling with porous septum. Maintenance-free. 1 m fixed cable with BNC+RCA plug
Description	Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature 0...60 °C
Diaphragm	Ceramic
Body diameter	12 mm
Electrolyte	Gel
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	5 m
Body material	Plastic
Name	201 T / 5
pH measuring range	0 ... 14 pH
temperature measurement range	0 ... +60 °C
Temperature sensor	Si / spina RCA Cinch
Reference system	Ag / AgCl