

Via della Meccanica, 25, 41012 Carpi (MO) - Italy
+39 059 653274 - info@giorgiobormac.com

VPT80-10H Cell Glass-Pt C=10+ATC range 100 to 500mS

Item Code: 50003372



DATA SHEET

Features

Cable	Si
Connection	BNC / RCA
	Platinum glass cell C=10, with built-in temperature sensor, fixed cable 1 m BNC-Cinch. Range 100µS...500 mS.
Description	ONLY for XS COND 80
Body diameter	12 mm
Image 2	https://app.onpage.it/api/storage/25ccad7f4b6ca27697e8cd44d317b6af3aed6
Body lenght	120 mm
Cable length	1 m
Body material	Glass
Name	VPT 10 H
Cond measuring range	100 µS/cm ... 500 mS/cm
temperature measurement range	0 ... +80 °C
Temperature sensor	Si / spina RCA Cinch